



**Options and Accessories Visuals**

**Base**

**Flat Free Tires - Block and V-Profile**



V-Profile for drive and front wheel for indoor use



Block-Profile for drive and front wheel for outdoor use

**V-Track Headrest (incl. Track extension)**



V-Trak Headrest winged  
 V-Trak Headrest simple

**V-Track Thoracic support & pad**



V-Trak Thoracic support & pad & cover large right & left  
 V-Trak Thoracic support & pad & cover medium right & left

**V-Track Backrest/Standard & Axxis**



Offers more back and neck support

**Control Options**

**R-Net Joystick Module**



R-Net includes connector points for accessory "Buddy Button"

**R-Net Omni+ Module**



Compatible with many input devices  
 Remote control  
 4/5 switch control  
 3 Axis switch or proportional head control  
 Sip & Puff  
 Single switch scanner  
 (For use with Flexible Arms or Joint Arm)

**Easy Touch**



Alternative keypress control with dual function positioning

**Swing away joystick holder**



Allows for easier transfer  
 Provides opportunity to maneuver chair in close proximity of tables

**Joystick protection guard**



Protects hands & Joystick

**LCEV Stand-up mechanism**



Standard on standing operation (armrest round with control switch)

## Control Options

### LCEV Armrest



Armrest 6 cm/2.4" wide with built-in control switch  
 Armrest 6 cm/2.4" wide with hand control  
 Armrest round with hand control  
 Armrest 5 cm/2" longer

### Swing away chest support



Rotatable and swing away to the rear. Increased feeling of safety for the user in standing position. Distance between backrest tube and the beginning radius

### LCEV Dual switch for table control



Stand-up operation with push buttons for table. Table tray not included

### Buddy Button



Green: ON/OFF  
 Yellow: MODE  
 Red: Emergency stop

### LCEV Control switch extension



Vertical for round armrest  
 Vertical for wide armrest  
 Horizontal

### Joystick-top



Customized control grips

T-bar    Mushroom  
 Stick    Chin  
 Softball    Ball

## Seating System

### Powered reclining back/leg rest 90-180°



Powers to complete prone/recline

### Powered recline and stand



Powers from lying position directly to standing

### Tilt in Space



Tilts 0-35 degrees

## Leg Rest Options

### Knee wedge

#### Rectangular



Knee supports are only delivered extra strong. Stabilization is removable

#### Conical



This knee wedges give more support, especially with spasticity

## Leg Rest Options

### Knee support - PRO



The new Knee support "PRO" option provides enhanced stability/rigidity and allows for separate positioning of each knee/leg

### Knee support - Integral



Individually adjustable

### One-piece footrest fix, no flip-up feature



Option- Easier to handle more stability

### Ground level footrest



Especially for users who use the wheelchair as a stand-up help

### Two-piece footrest



Standard- Easier to move in and out of chair

### Footplate protection frame



Protection frame only. Protects toes from collision injuries

### Footrest height adapter



Especially for users with short lower leg length. The footrest becomes 1.2" narrower than the seat width

### Seat depth / Footrest quick adjustment



This feature makes adjustment on footplate height quick and easy

Suitable for demo chairs

## Accessories

### Curb Climber



Front castor fork with double wheel. For curb climbing up to 3.2" / 8.1 cm

### Foam seat cushion with cover



Seat cushion option

### Backrest Extension - 6" (15CM)



Only available for standard Tension Adjust Sling Back backrest

### Back rest extension with push handle and upholstery



Length:  
5.9"  
9.8"  
13.8"

## Accessories

### Stroller handle/Push Handle Bar



Can be used with standard & V-Track back. Helps for manual pushing of the chair

### Table tray



Fixed tube in armrest

### Swing away table tray



Support installed at the armrest. Table is still fastened on the chair, even when not in use  
Only for Combi arm rest

### Swing away table tray center cut joystick holder



Depth adjustable. Table is still fastened on the chair, even if not in use. Joystick can be hinged away, if not used

### Skirt guard



Easy to remove for lateral transfer. (Pair)

### Thigh support



Supports the thigh and stabilizes the leg axle in seating and standing position

### Chest role



Individually adjustable in height. Offers additional support but does not replace the chest strap

### Foot guide



Can be adjusted individually. Secures and aligns feet

### Footrest wedge



Offers foot alignment for cases of spasticity and special need

### Chin control with shoulder holder



Shoulder- and chest holder with upholstery. Height & width are individually adjustable using "ON/OFF" & "MODE" button

### Chin control swing away SPZ



Electrically operated. Height, width & angle adjustable

### Attendant module



Attendant controls  
Easy operation transition

### Arms to adapt the R-Net Omni+ Module



Flexible arm including adapter  
Joint arm including adapter stronger and more adjustable  
Can be bent in every direction

### Monitor



Provides an objective outcome measurement for users and therapists  
Delivered with PC-Software.



## Accessories

### Rear wheels with drum breaks



Allows assistant to stop the wheelchair, safe and comfortably. Cannot be combined with the option of Spinergy, Camber or One arm drive

### Rear wheels "Spinergy"



Operation without additional energy. Black, blue, red, yellow & white

### Steel hand ring



Increased grip, less temperature-sensitive and less dirty hands

### Maxgrepp (Hand ring with rubber covered)



Original

Better hold for persons with limited strength or grip



Ergo Para CP

### One arm drive



For stroke patients or people with unilateral limited force or grip  
Only possible with standard wheel

### Push handle



Standard push handles

### Push handle extension adjustable



Helpful if the wheelchair is pushed by different persons

### Traditional Anti-tip rolls



Anti-tip rolls to keep chair from tipping backwards. Right & Left sold separately

### Anti-tip aid / tip aid swing away "HC"



Anti-tip aid to keep chair from tipping backwards which includes easy to use swing away function. Right & Left sold separately

### Brake lever extension (6 cm/2.4")



Simpler accessibility and less strength necessary

### Rear wheels with hub locking brakes



Light weight, improved. Cannot be combined with Hub locking brake, Drum brakes or One arm drive

### Spoke guard



Spoke guard-Transparent

Spoke guard-Transparent w/LEVO-Logo

## Vehicle Mount & Driving Accessories

### Vehicle tie-down kit



Secure tie-down for vehicle transport

### EZ-Lock brackets



Levo chairs adjust to EZ-Lock

### Lights and indicators (traffic approved)



Turn signals, brake lights, & lamps for on road use

### Rear view mirror



Mounted on Joystick side or mounted opposite side

### Mudguard for drive wheel



Mud Guards/Sold in pairs. Silver; Black & Rubber

### Frame reinforcement for e-fix



Wheelchairs with E-Motion need frame reinforcement. You can only have this combined with anti-tip aid "HC" and not with anti-tip rolls or tip aid