



5000™

easystand.com/5000



Renae Sommer
T6 Paraplegic, 5'5" (165 cm)

Who uses the 5000?

- Individuals from 5' to 6'5" (152-195 cm) and up to 280 lbs (114 kg)
- Individuals with Spinal Cord Injury, Multiple Sclerosis, Traumatic Brain Injury, Spina Bifida, Stroke & more
- Individuals with hip, knee or ankle contractures

What is the 5000?

The EasyStand 5000 is a supportive sit-to-stand stander designed to meet the needs of a variety of people. Tool-free adjustments and a wide variety of support options make it ideal for multi-user facilities.

Where would you use this stander?

- Rehabilitation Facilities and Hospitals
- Schools and Educational Settings
- Extended Care Units, Nursing Homes & Group Homes
- Home and Workplace
- Veterans' Hospitals



Rodney Kohout
L2 Paraplegic, 6'4" (192 cm)

When would you use this product?

- At home while reading or watching TV
- In rehab, when adjustability is key
- To get a good stretch
- Immediately after injury for the greatest benefits

"When I received the EasyStand I was standing within thirty minutes. It's a great product with simple assembly and ease of use – unassisted. My friends were amazed to see me at 6'2".

M.R. – Tampa, Florida, USA

"I've been a paraplegic now for 20 years & my EasyStand has become a very valuable tool in my life. Among the many things it has enabled me to do, it's allowed me to remodel 2 houses. How many paraplegics does it take to change a light bulb? One, with an EasyStand!"

S.W. – Sioux Falls, South Dakota, USA





How does it work?

- Transfer from the wheelchair to the seat.
- Position the kneepad and other supports.
- Pump actuator handle to desired standing position.

Why is the 5000 a better choice?

- Natural sit-to-stand transition allows user to stop at any position.
- User is fully supported from sitting to standing.
- It works well with people who have limited range of motion.
- Ease of use increases compliance with standing programs.
- Tool free adjustments and support options make the 5000 ideal for busy rehabs.



Steve Hansen
T5 Paraplegic, 6'0" (183 cm)

Adjustments

Quick & Easy, Tool-Free



Back Angle



Seat Depth



Table



Front Pad



Kneepad



Foot Plates
(plantar/dorsi, toe-in/toe-out)

"I like that the EasyStand is designed with the ability to add or remove a variety of options as needed."

J.N. –Tucson, Arizona, USA

"My son has a C6-7 injury – the EasyStand has made it possible for him to be able to move around the house while standing. It has also helped him strengthen his body."

T.M. – Spirit Lake, Iowa, USA

"There is a voice inside which speaks and says: This is the real me!"
William James



5000™ Options

easystand.com/5000-options




EasyStand® 5000™

The EasyStand 5000 comes standard with a padded seat & kneepad, foot plates, formed table with front pad, hydraulic actuator with handle, and graphite gray upholstery.

EasyStand® 5000XL™

The EasyStand 5000XL comes standard with all the features of the 5000. It is designed for people from 5' 5" (165 cm) up to 6' 10" (208 cm) and up to 350 lbs. (159 kg); seat depth range is 20"-26" (51-66 cm); seat to foot plate range is 13"-21" (33-53 cm); seat height is 21" (53 cm); Footprint is 31"x 37" (79x94 cm).

Upholstery Colors



#5200 Mobile Option
Provides independent mobility. Designed for the people who have good upper body strength. *Not available on 5000XL. Includes rear casters*



#5062 Locking Casters (set of four)
4" (10 cm) locking casters make it easy to roll the EasyStand from room to room on level surfaces.



#81019 Roho® Seat
Maximizes pressure relief. Features all the benefits of the Roho® Dry Floatation® technology.



#5020 Removable Actuator Handle
For people who tend to inappropriately lower themselves.



#5022 Actuator Handle Extension
For people with limited range of motion. Makes the first few pumps easier and releasing more convenient. *Adds additional 6" (15cm) to the reach*



#80634 Foot Straps (Pair)
Velcro® adjustment holds the feet in place. *Foot strap length over arch is 7"-17" (18-43cm)*



#80616 Secure Foot Option (Pair)
Velcro® and d-ring adjustment uses a four-point tightening system to hold feet in place. *Strap range over arch is 7"-17" (18-43 cm)*



#5026 Independent Kneepads
Accommodate knee contractures and leg length discrepancies; flip up for transferring. *Kneepad width center to center is 8" - 15" (20 - 38cm); kneepad depth range from front of seat is 3" - 8" (7-20 cm)*



#5040 Oversized Clear Table
Allows more room for reading and activities. *Not available with mobile option. Table size is 32"x 22" (81 x56 cm)*



#81011 Angle Adjustable Table (molded)
#80997 Angle Adjustable Table (oversized)
Attaches to existing table. Allows user to work at various angles from 10° to 45°.



#5030 Arm Troughs (Pair)
Adds support and positioning of the forearms for involved individuals; suction cups hold in place. *For use with the #5040 oversized clear table only*



#5250 No-table Option
Allows closer access to cabinets and other hard-to-reach areas. Works great with the mobile option.





#5024 Table Swing Away Sleeve
 Makes assisted transfers easier by moving the table away from the person.
Not available with mobile option



#5070 Ratchet Handles
 Can be used in place of existing knobs for ratchet motion and leverage.
Four small handles for table, two large for foot plates



#5038 Hip Guides (pair)
 Keep the body in alignment from sitting to standing. Works as a prompt to keep hips in alignment. *Width range is 11"-18" (28-46 cm); pad size is 5"x 6" (13x15 cm)*



#80246 Seat Angle Locator
 Helps therapists monitor progress. Attaches under seat with magnet.



#80809 Tool Pouch
 Attaches under the seat to hold the tools provided with the unit. *Tools not included*



#5100 Back Option
 For those who have limited upper body balance. *Height from seat is 11" (28 cm)*



#5100R Removable Back Option
 Allows the person to be transferred by a patient lift into the EasyStand. Also allows for more range of motion during therapy.
Back height from seat is 11" (28 cm)



#51500 Seat Belt
 Velcro® and d-ring adjustment provides support for user. *Fits hip circumference of 29"- 43" (74 - 109 cm)*



#5150 High Back Extension
 For people with minimal upper body balance. Includes rear anti-tips. *#5100 back option required. Height from seat is 22" (56 cm)*



#51600 Chest Strap - Large
 Velcro® and d-ring adjustment provides moderate upper body support. *Chest strap pad 7" W x 25"L (18W x 64L cm); Fits chest circumference of 35"-48" (89 - 122 cm)*



#80632 Chest Vest - Medium
#80633 Chest Vest - Large
 Neoprene® with camlock adjustment securely holds the person in position.
Medium size is 11"L x 9.5"W (28L x 24W cm); large size is 12.5"L x 10"W (32L x 26W cm)



#5036 Lateral Supports (pair)
 Aid in positioning a person who has limited trunk control. *Width range is 11"-18" (28 -46 cm); pad size is 5"x 6" (13 x15 cm)*



#5032 Head Support
 Works as a prompt to limit head movement.
Height range above seat is 20" - 31" (51-79 cm); depth range from back is 1" aft - 5" forward (3 cm aft - 13 cm forward)



#80643
 Universal Head Support Receiver
 Adapts to Whitmyer Biomechanix, Stealth® and Otto Bock™ headrests.

"The 5000 is the only model we've found that is easily adjustable and easy to use. It allows us to stand individuals according to their range of motion and ability. A wonderful and convenient product!"

M.D. – Berrien Springs, Michigan, USA

5000 Specifications

Height Range	5'0"-6'5"	152-195cm
Weight Limit	280 lbs.	114kg
Seat Depth Range (from seat pivot)	16"-24"	41-61cm
Seat to Foot Plate Range (from seat pivot)	11"-19"	28-48cm
Seat Height	19.5"	49cm
Table/Front Pad Depth Range (when standing)	2"-10"	5-25cm
Table Height Range (from seat pivot)	25.5"-35"	65-89cm
Kneepad Depth Range (from seat pivot)	3"-5"	8-13cm
Kneepad Width (center to center)	9"	23cm
Footprint	26.5"x33"	67x84cm
Weight of Basic Unit	118 lbs.	54kg
Frame Color	Black Powder Coat	
Standard Upholstery Color	Graphite Gray	
Black Molded Table Size	19.5"x20"	49x51cm

"Give yourself something to work toward-constantly."
 Mary Kay Ash