



Options and Accessories Visuals

7105 Northland Terrace Brooklyn Park, MN 55448 1-888-LEVO USA or 1-888-538-6872 Fax# 763-582-0442 amyj@danetechnolgies.com paulh@danetechnologies.com

Base

Flat Free Tires - Block and V-Profile



V-Profile for drive and front wheel for indoor use



V-Trak Headrest winged V-Trak Headrest simple



Block-Profile for drive and front wheel for outdoor use

V-Track Thoracic support & pad



Control Options R-Net Joystick Module

V-Trak Thoracic support & pad & cover large right & left V-Trak Thoracic support & pad & cover medium right & left

V-Track Backrest/Standard & Axxis



Offers more back and neck support

R-Net Omni+ Module



Easy Touch

R-Net includes connector points for accessory "Buddy Button"

Alternative keypress

function positioning

control with dual

The second

Compatible with many input devices Remote control 4/5 switch control 3 Axis switch or proportional head control Sip & Puff Single switch scanner (For use with Flexible Arms or Joint Arm)

Swing away joystick holder



Allows for easier transfer Provides opportunity to maneuver chair in close proximity of tables

Joystick protection guard



Protects hands & Joystick

LCEV Stand-up mechanism



Standard on standing operation (armrest round with control switch)

V-Track Headrest (incl. Track extension)

Control Options LCEV Armrest

Armrest 6 cm/2.4" wide with built-in control switch Armrest 6 cm/2.4" wide with hand

control Armrest round with hand control Armrest 5 cm/2" longer

LCEV Dual switch for table control



Stand-up operation with push buttons for table. Table tray not included

LCEV Control switch extension



Vertical for round armrest Vertical for wide armrest Horizontal



Swing away chest support

Buddy Button



Rotatable and swing away to the rear. Increased feeling of safety for the user in standing position. Distance between backrest tube and the beginning radius

> Green: ON/OFF Yellow: MODE Red: Emergency stop

Jovstick-top



Customized control grips

T-bar Mushroom Stick Chin Softball Ball

Seating System Powered reclining back/leg rest 90-180°



Powers to complete prone/recline



Powers from lying position directly to standing

Tilt in Space



Tilts 0-35 degrees

Leg Rest Options Knee wedge

Rectangular



Knee supports are only delivered extra strong. Stabilization is removable





This knee wedges give more support, especially with spasticity

Leg Rest Options Knee support - PRO



The new Knee support "PRO" option provides enhanced stability/rigidity and allows for separate positioning of each knee/leg

Knee support - Intregal



Individually adjustable

One-piece footrest fix, no flip-up feature



Option- Easier to handle more stability

Two-piece footrest



Standard- Easier to move in and out of chair

Footrest height adapter

Ground level footrest

Footplate protection frame



Seat depth / Footrest quick adjustment

Protection frame only. Protects toes from collision injuries

Especially for users who

use the wheelchair as a

stand-up help

Especially for users with short lower leg length. The footrest becomes 1.2" narrower than the seat width



This feature makes adjustmer on footplate height quick and easy

Suitable for demo chairs

Accessories Curb Climber



Front castor fork with double wheel. For curb climbing up to 3.2" / 8.1 cm

Backrest Extension - 6" (15CM)



Only available for standard Tension Adjust Sling Back backrest

Foam seat cushion with cover



Seat cushion option

Back rest extension with push handle and upholstery



Length: 5.9" 9.8" 13.8"

Accessories Stroller handle/Push Handle Bar



Swing away table tray



Support installed at the armrest. Table is still fastened on the chair, even when not in usé Only for Combi arm rest

Easy to remove for

lateral transfer. (Pair)

Can be used with standard & V-Track back. Helps for manual

pushing of the chair

cut joy stick holder

Swing away table tray center

Depth adjustable. Table is still fastened on the chair, even if not in use. Joystick can be hinged away, if not used

Fixed tube in armrest

Thigh support

Table tray



Supports the thigh and stabilizes the leg axle in seating and standing position

Chest role

Skirt guard



Individually adjustable in height. Offers additional support but does not replace the chest strap

Footrest wedge



Offers foot alignment for cases of spasticity and special need

Chin control swing away SPZ



Electrically operated. Height, width & angle adjustable

Arms to adapt the R-Net Omni+ Module



Flexible arm including adapter Joint arm including adapter stronger and more adjustable Can be bent in every direction



Monitor

Page 4 of 6



Provides an objective outcome measurement for users and therapists Delivered with PC-Software.

Can be adjusted individually. Secures and aligns feet

Chin control with shoulder holder



Shoulder- and chest holder with upholstery. Height & width are individually adjustable using "ON/OFF" & "MODE" button

Attendant controls

Easy operation transition

Attendant module

LEVO USA Options and Accessories Visuals 08/2010



Accessories Rear wheels with drum breaks



Steel hand ring

Allows assistant to stop the wheelchair, safe and comfortably. Cannot be combined with the option of Spinergy, Camber or One arm drive

Rear wheels "Spinergy"



Operation without additional energy. Black, blue, red, yellow & white

Maxgrepp (Hand ring with rubber covered)



Original

Better hold for persons with limited strength or grip

Standard push handles



Traditional Anti-tip rolls

H NONTO

Ergo Para CP

One arm drive



For stroke patients or people with unilateral limited force or grip Only possible with standard wheel

Increased grip, less temperaturesensitive and less dirty hands

Push handle extension adjustable



Helpful if the wheelchair is pushed by different persons

Anti-tip aid / tip aid swing away "HC"



Anti-tip aid to keep chair from tipping backwards which includes easy to use swing away function. Right & Left sold separately

Rear wheels with hub locking brakes



Light weight, improved. Cannot be combined with Hub locking brake, Drum brakes or One arm drive





Simpler accessibility and less strength necessary

Spoke guard



Spoke guard-Transparent

Spoke guard-Transparent w/LEVO-Logo

Anti-tip rolls to keep chair from tipping backwards. Right & Left sold separately

Vehicle Mount & Driving Accessories Vehicle tie-down kit



Secure tie-down for vehicle transport



Levo chairs adjust to EZ-Lock

Lights and indicators (traffic approved) Rear view mirror



Turn signals, brake lights, & lamps for on road use

Mudguard for drive wheel



Mud Guards/Sold in pairs. Silver; Black & Rubber



Mounted on Joystick side or mounted opposite side

Frame reinforcement for e-fix



Wheelchairs with E-Motion need frame reinforcement. You can only have this combined with anti-tip aid "HC" and not with anti-tip rolls or tip aid